



# QUICK *START* GUIDE





# CARE OF YOUR NEW TENNANT BIOMODULATOR®

***PLEASE KEEP THIS GUIDE HANDY  
AND REFER TO IT OFTEN***



# USING YOUR TENNANT BIOMODULATOR®



***PLEASE KEEP THIS GUIDE HANDY AND REFER TO IT AS NEEDED .  
YOU ARE ENCOURAGED TO ATTEND A COURSE WITH DR. TENNANT FOR  
MORE COMPLETE TRAINING. THE FOLLOWING ARE GENERAL SAFE USE  
GUIDELINES TO HELP YOU BEGIN ENJOYING BENEFITS.***



# Chapter 1 - Care of your new Tennant Biomodulator®

## *Table of Contents*

About this Guide.....	1
Statement of User Responsibility.....	2
Contraindications.....	3
Care of the Tennant Biomodulator®.....	4
Batteries.....	5
Powering on the Device.....	6
Accessory Port.....	7
Modes.....	8-12
Power Settings.....	13
How to Use Power.....	14
Colored Lights.....	15
Voltage Readings.....	16
Voltage Indications.....	17
Optional Attachments.....	18
Electrodes Attached with Wires.....	19-20
Warranty.....	21
FDA Classification.....	22



## Chapter 2 - Using your Tennant Biomodulator®

### *Table of Contents*

Powering on the Device.....	24
Practice Using Device.....	25-28
Treat Pain.....	29
Colored Lights.....	30
Power Settings.....	31
How to Use Power.....	32
Voltage Readings.....	33
Modes.....	34
Voltage Indications.....	35
Measure BioTerminals.....	36-39
Balancing Autonomic Nervous System.....	40-43
Using Biomodulator for Extended Time & Sleep Mode.....	44
Electrodes Attached with Wires.....	45
Improve Total Body Voltage.....	46-47
Connecting Black and Red Leads.....	48
Troubleshooting Voltage Readings.....	49
Simple Daily Treatment.....	50
Things to Remember.....	51
Register for Health Conference.....	52





# About this Users Guide

- Please read all information in this guide and refer to it often.
- This material is only intended to provide initial instructions about the basics of turning your device on and understanding its features / functionality.
- Please carefully read the **Contraindications** on page 3 before using the device on anyone.
- You are encouraged to attend a formal training course taught by a licensed medical professional, so you will gain the full potential for achieving desired results as soon as possible.
- You should **NOT** use your device in a professional therapy setting without having been properly trained.
- Technical support offered by the distributor is for questions about the mechanics of the device.
- Technical support does **NOT** include instruction about the medical use of the device nor to give medical advice.

# Statement of User Responsibility

- You have obtained your **Tennant Biomodulator®** either because you are a licensed healthcare professional or because you have a prescription for the device on file with Senergy Medical Group.
- The use and care of the device is solely your responsibility.
- Examples provided by other users are for informational purposes only and are not a substitute for the advice provided by your own physician or any other medical professional.
- Examples are to give you more information with which to make informed decisions.



# Contraindications

- Do not use with **pacemaker** or **implanted electronic device** without permission of cardiologist.
- Do not use on a **stroke patient** until a possible blood clot has had time to resolve.
- Do not rub device over a vessel that might contain a **blood clot**, e.g., a **tender varicose vein** or the **carotid arteries**.
- Do not use when someone is under the **influence of alcohol or street drugs**.
- Do not use electrode directly on **eyeball**.
- Do not use **during pregnancy** (for liability reasons)
- Do not use on a patient who is very **sensitive to electromagnetic frequencies**, e.g., someone who can't be around a computer.
- Do not use on a patient who **doesn't want you to use it**.

# Care of Device

Your Tennant Biomodulator® has software and advanced circuitry like a computer, it is important to remember this when handling your device.

- Don't turn on if below freezing.
- Clean electrode only with a 70% isopropyl alcohol .
- Do not leave your device in a hot car in the summertime.
- **Do not try to open your device; doing so will void your warranty.**
- Device can be damaged to the point of non-function if dropped on hard surface. (a protective carrying case is recommended)
- Device WILL be damaged beyond repair if dropped in water or liquids are spilled inside.
- The battery compartment slides off and on the device. Do not force the battery cover.
- Be careful when accessories are plugged into the port. Damage can occur inside the port if there is too much movement at the plug.
- Always unplug wires from the accessory port and store device in protective case.



Electrode corrosion can occur  
If not cleaned regularly.



# Batteries

- Tennant Biomodulator® uses 2 AA alkaline batteries included with the device.
- Do Not mix old and new batteries, as this can cause them to leak.
- When ON/OFF switch is in the ON position, device is still powered and battery life is being used.
- The device is programmed to go into sleep mode after one hour, to save battery life. To wake it up again push the (+) button.
- Device will signal you that the batteries are running low by a series of consecutive beeps or by switching back to the Assess mode.
- Please do not call us to say that your device is broken because it won't stop beeping.
- If you replaced the batteries and it is still beeping, go buy a new package of batteries and replace both batteries with new ones.
- If you return your device and it is determined that it just needed new batteries you will be charged \$50.

# Turn Tennant Biomodulator® On

The **ON-OFF** switch is on the side of the unit. Slide the switch up to turn it on. Slide the switch down to turn it off.




# Optional Accessory Port

Optional **accessory port** on top right corner of unit.

This port should be used for authorized accessories only!

Do NOT force anything into this port or the warranty may be voided.

**Always unplug wire when not in use.**



# Mode Selection with Tennant Biomodulator®

The illuminated **BLUE** light indicates the active mode.



Each time you press the **MODE** button, the illuminated blue light will move to the selected mode.



# Assess Mode

- **Assess** – “voltmeter” Reads the voltage when placed on the skin. This is used to measure specific points on the body. Place the metal contact of the device on the skin to be measured. A single light will illuminate to indicate if voltage is **normal, healing or low**. Device must remain still. Once the light is visible, device can be moved to assess another area. (convenient when you want to record your voltage on a chart)

# Ten-8™ mode

- **Ten-8™** – A therapy mode specifically developed for relief of chronic, intractable pain and for muscle and bone. Device does not have to be stationary; can glide over areas of concern.

**This mode does not read voltage.** The sensation felt in Ten-8™ is consistent in intensity. The colored lights are indicating power level.

# Infinity™ mode

- **Infinity™** - A treatment mode designed for symptomatic relief of pain and for organs. Device does not have to be stationary; can glide over the area of concern. **This mode does not read voltage.** The sensation felt in Infinity™ starts out low then it increases. It goes from low to high in intensity. The colored lights are indicating power level.

# Automatic mode

- Automatic – This is a two-phase mode: (1) Reads your voltage and then (2) Infinity™ therapy mode begins automatically and continues 1 minute. A series of beeps sound when therapy session ends. It will then take another reading, so you can see if there was any improvement. Either move and treat another location or remain stationary for the next treatment cycle.
- Use Automatic when you want to measure the voltage and immediately treat.
- The goal is to bring up your initial reaction to above 25. This should happen after 2 or 3 minutes.
- If initial reaction is not improved in 2 – 3 minutes put the device in a treatment mode and treat for several minutes.
- **Not intended to be used with the electrode pads.**
- **Not intended to be used for treatment longer than a few minutes because it resets every minute.**

# Tennant Biomodulator® Power Settings

## when in the Ten-8™ or Infinity™ modes

### Power Setting Indicators

#### When in Ten-8™ or Infinity™

The lights on the curve will illuminate as the power of the unit is increased.

As the “+” button is pressed and held, the lights will start to flash, then go solid, and the next light in line will begin to flash as the power increases.

The lights will continue to illuminate in order, from left to right.

Once at maximum power, the last light on the curve will be solid orange (10 power).



### Adjusting the Power Levels



Press & hold the **UP** (+) button to Increase

Power (push and hold down the button to control power)



Press & hold the **DOWN** (-) button to decrease

power (push or hold down the button to control power)

**NOTE – Device indicates power level when in Ten-8™ or Infinity™ modes only. When in Assess or Automatic the device is indicating a voltage reading.**

# Setting Power Level

- You set the power to what feels comfortable
- A slight tingle should be felt
- It is better to use **too little power** (rather than too much).

**The cells will shut down if you use too much power.**

**(This is not harmful but it is counterproductive.)**

- Let children control the power themselves. Be sure to use low power for animals.
- When in **Ten-8™ or Infinity™ modes**, the illuminated lights ranging from green to orange are indicating only the power settings, and are **NOT reading voltage.**
- When in **Ten-8™ mode** **the sensation felt will be consistent in intensity.**
- When in **Infinity™ mode** **the sensation felt starts out low then it increases. It goes from low to high in intensity. (use caution and increase power slowly)**

**Notice the colored lights indicate something different if you are in an assessment mode or a therapy mode.**

### *Assessment Modes*

- If you are in Assess or Automatic modes the light is indicating a voltage reading.

### *Therapy Modes*

- If you are in Ten-8™ or Infinity™ modes the light is indicating the power level.



# Voltage Readings in Assess or Automatic Mode

*Voltage readings can only be taken in Assess or Automatic*

- Use Assess when you want to record your voltage on a chart.
- Use Automatic when you want to measure the voltage and immediately treat.



*When in Assess or Automatic Modes the light is indicating a voltage reading not power level.*



# Voltage Indications

## When in Assess or Automatic Modes



# Optional Attachments

There is an Optional Accessory Port on the other side of the unit.

There are many different attachments that can be plugged into the

**Tennant Biomodulator®** that make it much easier to treat specific areas of the body.



- **Scalp Electrode**

- *For use in hair on scalp, chest, etc.*
- *Can also be used on animals*



- **Face Electrode**

- *Necessary for taking readings of voltage on ones own face or back*
- *For general use in small areas*



- **Pencil Electrode**

- *For use on acupressure points*
- *Used inside the mouth for problems with teeth*



- **Large Body Electrode**

- *Has an 8" handle, great to help reach areas like the back*
- *Useful for practitioners too*

Replaces direct use of device by extending reach of electrode.

The **Tennant Biomodulator®** is a stand alone device these attachments are optional.

**Always unplug any attachment when not in use.**

# Electrodes Attached with Wires



The electrode pads, can be placed anywhere on the body, whenever an extended therapy is needed.



The electrode sleeve and glove are helpful when a therapy is needed on rounded or bony surfaces.



# Electrodes Attached with Wires



The electrode pads, and electrode sleeve, sock or glove are connected with a wire. The wire splits into two parts one side has a red end and the other side a black end that plugs into these items. When using the electrode pads you would use 2 pads, or 4 depending on which wire you use.



*Use this wire to connect 2 pads or a sock or glove*



*Use this wire to connect 4 pads at once*

# One-Year Limited Manufacturer Warranty

- The warranty begins on the date of purchase and covers manufacturer defects only.
- The warranty becomes invalid upon any of the following:
  - More than one year has passed since the purchase date.
  - The device was damaged or broken due to abuse, misuse or improper use.
    - Dropping it, throwing it, leaving it in extreme heat or cold, etc.
  - Water or chemical damage.
  - Broken seal or case; trying to open device.
  - Shorting of electrodes.
  - Violent vibration.
    - NOTE: violent vibration may damage components and reduce the effectiveness of the device.

# FDA Classification

## Important Notice

FDA Classification: 882.5890 Neurology transcutaneous electrical nerve stimulator for pain relief

---

882.5850 Neurology biofeedback device

---

Identification for Use: Symptomatic relief and management of chronic, intractable pain

---

Adjunctive treatment in the management of post-traumatic surgical and post-traumatic pain

---

Relaxation training and muscle reeducation

---

Caution:



# **In the next chapter we will be helping you to get started using your new Tennant Biomodulator®**

- We want you to feel comfortable using your device immediately.
- Dr. Tennant encourages you to read this material watch the DVDs and start using your Biomodulator right away.
- Dr. Tennant feels you will get much more from his training course, if you have been using it for a while.
- In this next chapter we are giving you a head start on some things that Dr. Tennant teaches in the course.





# Turn Tennant Biomodulator® On

The **ON-OFF** switch is on the side of the unit. Slide the switch up to turn it on. Slide the switch down to turn it off.



# Put The Device in Ten-8™ Mode By Pressing the Mode Button

The illuminated **BLUE** light indicates the active mode.



Each time you press the **MODE** button, the illuminated blue light will move to the selected mode.

# Put The Device On Your Skin And Press The + Button Until You Feel A Slight Tingle


- Move the device around on your skin to get an understanding of what it feels like.
- It is important to feel comfortable with using your Biomodulator.
- Glide it all over the skin.
- Glide it up and down your arms, legs, torso, hands, feet.



Using your device in this way provides a general benefit.

# Do Not Be Overly Concerned About Making Mistakes

- Some users feel initial concern about making mistakes. When you follow the instructions you need not be concerned.
- The Biomodulator was designed for general use and does not have to be used precisely.
- Even if you are providing the therapy to a different area than intended, it is not harmful and, in most cases, will still benefit you.
- There are many good ways to accomplish delivery of additional energy – do not be overly concerned about making mistakes.



# Practice Using Your Device In Both The Ten-8™ Mode And The Infinity™ Mode

- Feel the difference between the two modes.
- It is important to understand that there are guidelines for determining which mode may be optimal for different situations. However there is no harm in experimenting.
- **These guidelines are -**
- **Ten-8™** is used for muscle and bone.
- **Infinity™** for all other issues related to organs and general health.





# Treat Pain As Long As Needed

- Treat any areas of pain by placing the device on the painful area.
- The device can be held on one spot or it can be rubbed over the area.
- The device can be used as long as needed.
- If the pain is from an injury of bone or muscle **Ten-8™** is the mode of choice.
- If the pain is from a stomach ach **Infinity™** would be the mode of choice because it is an organ.



# Notice The Colored Lights Indicate Something Different If You Are In An Assessment Mode or a Therapy Mode

## Assessment Modes

- If you are in Assess or Automatic modes the light is indicating a voltage reading.

## Therapy Modes

- If you are in Ten-8™ or Infinity™ modes the light is indicating the power level.



# Tennant Biomodulator® Power Settings

## when in the Ten-8™ or Infinity™ modes

### Power Setting Indicators

#### When in Ten-8™ or Infinity™

The lights on the curve will illuminate as the power of the unit is increased.

As the “+” button is depressed and held, the lights will start to flash, then go solid, and the next light in line will begin to flash as the power increases.

The lights will continue to illuminate in order, from left to right.

Once at maximum power, the last light on the curve will be solid orange (10 power).



### Adjusting the Power Levels



Depress the **UP** (+) button to Increase Power (push and hold down the button to control power)



Depress the **DOWN** (-) button to decrease power (push or hold down the button to control power)

**NOTE – Device indicates power level when in Ten-8™ or Infinity™ modes only. When in Assess or Automatic the device is indicating a voltage reading.**



# How to Use Power

- You set the power to what feels comfortable
- A slight tingle should be felt
- It is better to use **too little power** (rather than too much).

**The cells will shut down if you use too much power.**

**(This is not harmful but it is counterproductive.)**

- Let children control the power themselves. Be sure to use low power for animals.
- When in **Ten-8™ or Infinity™ modes**, the illuminated lights ranging from green to orange are indicating only the power settings, and are **NOT reading voltage.**
- When in **Ten-8™ mode** **the sensation felt will be consistent in intensity.**
- When in **Infinity™ mode** **the sensation felt starts out low then it increases. It goes from low to high in intensity. (use caution and increase power slowly)**

# Voltage Readings in Assess or Automatic Mode

*Voltage readings can only be taken in Assess or Automatic*

- Use Assess when you want to record your voltage on a chart.
- Use Automatic when you want to measure the voltage and immediately treat.



*When in Assess or Automatic Modes the light is indicating a voltage reading not power level.*

# Modes

- **Assess** – “voltmeter” Reads the voltage when placed on the skin. This is used to measure specific points on the body. Place the metal contact of the device on the skin to be measured. A single light will illuminate to indicate if voltage is **normal, healing or low**. Device must remain still. Once the light is visible, device can be moved to assess another area. (convenient when you want to record your voltage on a chart)
- **Ten-8™** – A treatment mode specifically developed for relief of chronic, intractable pain and for muscle and bone. Device does not have to be stationary; can glide over areas of concern. **This mode does not read voltage.** **The sensation felt in Ten-8™ is consistent in intensity. The colored lights are indicating power level.**
- **Infinity™** - A treatment mode designed for symptomatic relief of pain and for organs. Device does not have to be stationary; can glide over the area of concern. **This mode does not read voltage.** **The sensation felt in Infinity™ starts out low then it increases. It goes from low to high in intensity. The colored lights are indicating power level.**
- **Automatic** – This is a two-phase mode: (1) Reads your voltage and then (2) **Infinity™** therapy mode begins automatically and continues 1 minute. A series of beeps sound when therapy session ends. It will then take another reading, so you can see if there was any improvement. Either move and treat another location or remain stationary for the next treatment cycle. Use Automatic when you want to measure the voltage and immediately treat.



# Voltage Indications

## When in Assess or Automatic Modes





*When in Assess or Automatic Modes, the light is indicating a voltage reading not power level. Each light indicates voltage range.*

# Measuring and Treating BioTerminal Points

- When you want to take voltage readings place your device firmly on each point shown in the image on the next page.
- Use the **Assess mode** if you want to write down the readings and come back to do treatment later. This is helpful if you want to chart your results.
- Treat each point for 3 minutes in the **Infinity™ mode**.
- If you want to quickly **measure voltage & treat**; use the **Automatic mode**.
- Automatic takes a reading and then automatically starts a one minute treatment in the **Infinity™ mode**. You will hear a series of beeps and then it will take another reading. It will continue this cycle of taking readings and giving a one minute treatment.
- It is not necessary to take measurements every day.

# Measure Voltage on each BioTerminal Point



You can make copies of this chart for future use

# Avoid Measuring Voltage When

- The skin is hot and sweaty
- The skin is very dry
- You have not had at least 8oz. of water to drink in the last 30 min.

If you take readings under these conditions, they will not be accurate

- **If the skin is too dry the reading will be lower**
- **If the skin is sweaty the reading will be higher**

# Only Measure Voltage Once A Day

- Skin should be clean and dry free from makeup or moisturizers.
- When measuring your points, it is important to place the electrode firmly on that point.
- The first reading you get is the one that counts.
- Please note that your voltage can be different from day to day.
- When we exercise our voltage will be higher because our muscles act like rechargeable batteries and hold on to the voltage.
- It is very important to drink at least 8<sub>oz.</sub> of water 30 minutes before measuring voltage.
- Taking measurements is to help you gather information it is **not something that you have to do every day.**



# **“Tennant Autonomic Balancing System” Balancing the Autonomic Nervous System**

**“Fixing Your Bowling Ball”**

**Dr. Tennant teaches the importance of correcting imbalances in the body.**

**The effects that take place with balancing of the autonomic nervous system are both physically seen and some unseen, but most definitely felt.**

**This treatment is very important because when your autonomic nervous system is out of balance it effects everything, so we do this technique every day to stay in balance.**

**This is a key factor in chronic pain and chronic illness.**

Examples of physical effects that may be noticed:

- The two sides of the face may look asymmetrical.
- The jaw can move to one side, causing TMJ.
- Sinus and Nasal obstruction which can lead to snoring.
- Vertebrae C1 and C2 move to one side causing headaches and neck aches as well as migraine headaches.
- The spine can curve causing extrusion of disks.
- One shoulder may be higher than the other, making one arm seem short.
- The pelvis can rotate giving low back pain and disk extrusion and possibly making one leg shorter than the other. This places more weight on one hip-knee-ankle making those joints wear out.
- The locking of the cranio-sacral pump causes the entire nervous system to use stagnant cerebrospinal fluid resulting in a general decrease in its function.

# “Tennant Autonomic Balancing System” Balancing the Autonomic Nervous System


“Fixing Your Bowling Ball”

**This is done by using the Tennant Biomodulator® in the Infinity mode.  
( Daily Treatment )**

1. Place the Biomodulator on the side of the neck half-way between the shoulder and the ear. (See image on next page.) Unit should be in **Infinity™ mode** and power should be set at a comfortable level. Be prepared to increase or decrease power as treatment continues.
2. You will need to push firmly into this muscle, as you are trying to reach the nerve behind the muscle.
3. When you find the right spot you will know, because it feels different. You should feel a tingling sensation along the back of the ear. It may take time to find the right spot so be patient.
4. Hold unit in place for 2 – 3 minutes.
5. Repeat treatment on other side of neck.
6. Once both sides are complete, turn head left to right, approximately 10 times each side.
7. Rotate shoulders front and back, approximately 10 times each way.
8. Stand and rotate at the waist, approximately 10 times.
9. Walk around room or area for a short distance.

# This Image Shows Where to Place the Biomodulator for Balancing the Autonomic Nervous System

*“Tennant Autonomic Balancing System”  
“Fixing Your Bowling Ball”*



# “Tennant Autonomic Balancing System”

## Balancing the Autonomic Nervous System

We want to further explain the bowling ball technique because it is so important. This is a learned technique and it will take time to get it just right.

**Balancing the autonomic nervous system is something that will need to be done every day because it is important to keep the autonomic nervous system in balance.**

In the video Dr. Tennant is demonstrating the technique and saying to put the device in the Automatic mode because it will beep after just one minute to let you know that you are finished.

We have found that it is easier for people who are just learning to do this to use the Infinity mode instead. The reason is we don't want you to limit yourself to just one minute. It may take you one minute to find the right spot. When you find the right spot you will know because it feels differently. You should start to feel the sensation radiating in different areas. It feels very strange, a tingling sensation along the back of the ear. Once you feel this strange sensation you know you are doing it right. Now that you have felt that tingle in the ear keep the Biomodulator in place for at least one minute. You can do this for several minutes if you want to. The longer you do this the more relaxed you become.

In the video some have misunderstood and thought that Dr. Tennant is saying that you should not feel anything. That is not correct he is simply saying that the Infinity mode modulates higher in intensity and not to have the power too high. **This technique is the most important thing for you to learn to do correctly and continue to do every day.**

# Using Your Tennant Biomodulator® for Extended Time.

- You may want to use your Biomodulator for an extended amount of time for things like back pain and knee pain.
- The electrode pads are very convenient because they stick to the body and can be worn under clothing.
- The electrode pads can be used as long as is needed. (there is no time limit)
- **When wanting to use your device more than one hour you need to bypass the sleep mode**
- **Bypass Sleep mode** – To temporarily bypass the sleep mode  
Before turning on device, you must hold down both (+) and (-) buttons at the same time, then turn on the device while still holding down the (+) and (-) buttons. Count to 3 before you release the (+) and (-) buttons. This bypass of the sleep program is only temporary. When you turn the device off again it will go back to the default programming. (This will bypass the sleep mode & silence the device)
- **Make sure you are not using the Automatic Mode with the electrode pads.**

# Electrodes Attached With Wires



The electrode pads, and electrode sleeve, sock or glove are connected with a wire. The wire splits into two parts one side has a red end and the other side a black end that plugs into these items. When using the electrode pads you would use 2 pads, or 4 depending on which wire you use.



*Use this wire to connect 2 pads or a sock or glove*



*Use this wire to connect 4 pads at once*

# Improve Total Body Voltage

Dr. Tennant teaches how to improve total body voltage.

- This is done by placing the electrode pads at specific points on the body and keeping them on for many hours. This is recommended for people with a chronic illness or just low energy.
- You would use the **Infinity™** mode and keep the pads in place as long as needed.
- The image on the next page is showing how you can improve your body's total voltage by placing electrode pads at each of these points.

**If you want to connect all 4 locations at once,  
you will need to use the connection with 2 black and 2 red ends.**

**\*This can be done all night while you sleep.\***

***Remember to set your Biomodulator to stay on all night and keep it quiet  
you must hold down the + and - buttons while turning on the device.***



# This Image Shows Where to Place the Electrode Pads to Improve Total Body Voltage.

- Dr. Tennant teaches this placement to help a person achieve benefit throughout their body.
- Attach the black & red wires as shown in the picture.
- Black on the front and red on the back.
- When working on the **midline** of the body, as is shown in this image, always put the black lead on the front and the red lead on the back.






# Place Black Lead Wherever You See “Black”

# Place Red Lead Wherever You See “Red”

If you are using hand grips hold the red in the left hand and the black in the right hand.

- This image is showing where to place the electrode pads to keep optimal voltage flowing, when using pads on different quadrants of the body, arms, legs etc.
- The pads can be placed left to right, or front to back.
- Dr. Tennant teaches this placement because of how voltage is in our bodies. Following this guide should get you good results.



Remember when working on the midline always put the black in front and the red in back as shown on previous image.

# Troubleshooting Voltage Readings

- If you can not get a reading; resolve these issues first
  - Skin is too dry – drink more water, take readings about 30 minutes after a bath or shower when skin will have more moisture. But do not add water or lotion to the skin in order to moisten skin.
  - Person is dehydrated – drink more water
  - Reading takes longer to display (may take up to 10 seconds) when voltage is very low
  - At times the voltage can be too low to be able to get a reading because of a serious condition – start a treatment by using **Ten-8™** and **Infinity™** modes. You will notice when voltage is especially low you will need to increase the power, and as voltage is improved in the area, you will need less power, (you continue to adjust the power as needed to keep it at a comfortable level)

# Simple Daily Treatment

- Drink 8<sub>oz</sub> of pure alkaline water
- Balance your autonomic nervous system “fixing your bowling ball”
- Give your body energy by either putting the device on each BioTerminal point in the Infinity mode for 3 min. each or using the hand grips for about 10 min.
- Now that you have balanced your body and provided a base of energy for your body you can work on any areas of pain or chronic illness.
- An easy way to work on chronic pain or illness is with the BioTransducer, simply point the light at any area of concern.
- Another way to work on chronic pain or illness is to use the electrode pads. You can put them on and wear them all day if you want. See images on pages 47 & 48.

# Things to Remember

- **Infinity™ mode** – starts out low and increases in intensity, wait and increase the power slowly and adjust for comfort.
- **Sleep mode** - Device is programmed to go into sleep mode after one hour of use; to wake it up again push the (+) button.
- **Bypass Sleep mode** – To temporarily bypass the sleep mode

Before turning on device, you must hold down both (+) and (-) buttons at the same time, then turn on the device while still holding down the (+) and (-) buttons. Count to 3 before you release the (+) and (-) buttons. This bypass of the sleep program is only temporary. When you turn the device off again it will go back to the default programming. (This will bypass the sleep mode & silence the device)

# Register for Dr. Tennant's Health Conference & Biomodulator Training

- If you are not yet registered to attend a conference, please contact Senergy Medical Group to do so.
- We value your comments and appreciate your sharing these with us.
- Hope to see you at a conference soon!



Senergy Medical Group

9901 Valley Ranch Parkway East · Suite 1009 · Irving, TX 75063

972-580-0545 · 866-514-8221 Toll Free · 214-260-7765 Fax

[support@senergy.us](mailto:support@senergy.us) · [www.senergy.us](http://www.senergy.us)